
# 14. Cromatografia

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### Basic principles of SPE


Solid phase extraction (SPE) is a powerful method for sample preparation and is used by most chromatographers today.

More than 20 years ago MACHEREY-NAGEL designed and introduced CHROMABOND® SPE cartridges containing silica-based adsorbents. Since then we developed the widest range of phases and products for SPE based on silica and polymeric materials.

SPE has capabilities in a broad range of applications:

- environmental analyses
- pharmaceutical and biochemical analyses
- organic chemistry
- food analysis

## Solid phase extraction


SPE is a form of digital (step-wise) chromatography designed to extract, partition, and/or adsorb one or more components from a liquid phase (sample) onto a stationary phase (adsorbent or resin). An adsorbed substance can be removed from the adsorbent by step-wise increase of elution strength of the eluent (step gradient technique). SPE extends a chromatographic system's lifetime, improves qualitative and quantitative analysis, and the demand placed on an analytical instrument is considerably lessened.

#### In general, SPE is used for three important purposes in state-of-the-art analyses:

- concentration of the analyte (up to factor 10.000 - increase of chromatographic sensibility / improved limits of detection)
- removal of interfering compounds (protection of subsequent analyses like HPLC, GC, TLC, UV or IR spectroscopy, ...)
- changing an analyte's environment to a simpler matrix more suitable for subsequent analyses

#### Advantages of SPE compared to classical liquid-liquid extraction:

- lower consumption of solvents
- faster - enormous time savings
- lower costs per sample
- potential for automation
- high consistency in individual sample handling
- more specific selectivity because of the broad range of adsorbents and different retention mechanisms
- optimisation of extraction by variation or adjusting of the solid phase and chromatographic conditions


### Basic principles and summary of MN phases for SPE

#### Design of columns and cartridges

All CHROMABOND® columns and CHROMADIX® cartridges are manufactured from polypropylene (PP) with lowest content of extractables (plasticizers, stabilisers, ...) offering blank value free results by usage of most common solvents. The high quality CHROMABOND® adsorbents are kept in place by chemically very inert polyethylene filter elements (PE, standard pore size 20 µm). Funnel-shaped large volume columns, 96-well plates and cartridges for automated and on-line SPE packed with CHROMABOND® adsorbents are available on request.

#### CHROMABOND® polypropylene columns

- PP columns with PE filter elements
- different sizes from 1, 3, 6 up to 150 ml
- adsorbent weights from 20 mg to 50 g
- male luer tip as exit
- compatible with most robots (e.g., Gilson ASPEC™, Caliper AutoTrace®)

CHROMABOND® glass columns are available on request.

#### CHROMAFIX® cartridges

- PP cartridges with PE filter elements
- three different sizes with different adsorbent weights: **Small** (0.4 ml), **Medium** (0.8 ml), **Large** (1.8 ml)
- female Luer tip at the inlet, male Luer tip as exit
- offers alternative way of handling using positive pressure by syringes or peristaltic pumps
- especially suited for convenient SPE of small sample volumes

## Solid phase extraction

Phase	Matrix	Modification / Application	Similar phases*
<b>RP phases</b>			
HR-X	PS/DVB		ENVI-Chrom P, Strata™-X, Oasis® HLB, Nexus
PS-RP	PS/DVB	removal of organic components	Strata™ SDB-L, Bond Elut® ENV, Bond Elut® LMS, DCS-PS/DVB, ENV PS-DVB, Bakerbond™ H <sub>2</sub> O-phobic DVB, Isolute® 101
C <sub>18</sub>	silica	octadecyl, not endcapped	Strata™ C18-U, Accubond® C18, Bakerbond™ PolarPlus, Isolute® C18, LiChrolut® RP-18
C <sub>18</sub> Hydra	silica	octadecyl, not endcapped, for polar analytes	
<b>Normal phases</b>			
SiOH	silica	unmodified	Strata™ Si-1, Bond Elut® silica, DSC-Si, LC-Si, CLEAN-UP® silica, Accubond® silica, Bakerbond™ silica gel, Isolute® silica, LiChrolut® Si
NH <sub>2</sub>	silica	aminopropyl	Strata™ NH <sub>2</sub> , Sep-Pak® NH <sub>2</sub> , Bond Elut NH <sub>2</sub> , DSC-NH <sub>2</sub> , LC-NH <sub>2</sub> , CLEAN-UP® aminopropyl, Accubond® NH <sub>2</sub> , Bakerbond™ amino, Isolute® NH <sub>2</sub> , LiChrolut® NH <sub>2</sub>
OH	silica	diol	DSC-Diol, LC-Diol, Accubond® Diol (OH)
Alox A/N/B	aluminium oxide	acidic / neutral / basic	LC-Alumina-A/N/B, Accubond® aluminium oxide A/N/B
Florisil®	magnesium silicate		Strata™ FL-PR, Sep-Pak® Florisil®, Bond Elut® Florisil®, ENVI-Florisil®, LC-Florisil®, CLEAN-UP® Florisil®, Accubond® Florisil®, Bakerbond™ Florisil®, Isolute® FL, LiChrolut® Florisil®
PA	polyamide 6		DPA-6S
<b>Ion exchangers</b>			
SA	silica	benzenesulphonic acid cation exchanger (SCX)	Strata™ SCX, Bond Elut® SCX, DSC-SCX, LC-SCX, CLEAN-UP® Benzenesulfonic Acid, Accubond® SCX, Bakerbond™ Aromatic Sulfonic Acid, Isolute® SCX, LiChrolut® SCX
PSA	silica	propylsulphonic acid cation exchanger	
PS-OH <sup>-</sup>	PS/DVB	strong anion exchanger, OH <sup>-</sup> form	Oasis® MAX
PS-H <sup>+</sup>	PS/DVB	strong cation exchanger, H <sup>+</sup> form	Oasis® MCX, Strata™ X-C
PS-Ag <sup>+</sup>	PS/DVB	strong cation exchanger, Ag <sup>+</sup> form	
PS-Ba <sup>2+</sup>	PS/DVB	strong cation exchanger, Ba <sup>2+</sup> form	
<b>Phases for special applications</b>			
Drug	silica	bifunctional C <sub>8</sub> /SA, for enrichment of drugs from urine	Strata™ Screen-C, Bond Elut® Certify I, DSC-MCAX, Clean Screen® DAU, Accubond® Evidex, Bakerbond™ Narc-2, Isolute® HXC, LiChrolut® TSC
Drug II	silica	bifunctional C <sub>8</sub> /SB, for extraction of THC and derivatives as well as acidic analytes from biological fluids	Strata™ Screen-A, Bond Elut Certify II, Clean Screen® THC, Bakerbond™ Narc-1, Isolute® HAX
CN/SiOH	silica	combination phase for PAHs from soil	
NH <sub>2</sub> /C <sub>18</sub>	silica	combination phase for PAHs from water	
Na <sub>2</sub> SO <sub>4</sub> /Florisil®	silica	combination phase for extraction of hydrocarbons from water (DIN H-53 / ISO DIS 9377-4)	
SA/SiOH	silica	combination phase for PCBs from waste oil	
Diamino	silica	primary and secondary amine functions (PSA), for de-termination of pesticides in food (QuEChERS method)	Supelclean PSA, Bond Elut PSA

\* phases which provide a similar selectivity based on chemical or physical properties (list not complete)

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
#### Standard protocol for method development with CHROMABOND® HR-X

Column type: CHROMABOND® HR-X  
3 ml / 200 mg

Sample pretreatment: if necessary, adjust pH value

Column conditioning: 5 ml methanol

Equilibration: 5 ml water

Sample application: slowly aspirate the sample through the column

Column washing: 5 ml water - methanol (95:5, v/v)


Elution: after drying 3 x 2 ml methanol

Further analysis: if necessary, evaporate and redissolve in a suitable solvent; HPLC or GC

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
2

### PAHs from water containing humic acids

**Column type:**  
CHROMABOND® NH<sub>2</sub>/C<sub>18</sub>, 6 ml, 500 mg/1 g  
glass column

**Sample pretreatment:**

mix 500 ml water sample with 25 ml 2-propanol

**Column conditioning:** 10 ml dichloromethane, 10 ml methanol, then 10 ml dist. water – 2-propanol (9:1, v/v)

**Sample application:** aspirate 500 ml of the pretreated water sample through the column (~ 5 ml/min)

**Washing:** 2 ml dist. water – 2-propanol (9:1, v/v), then dry column (about 20 min, vacuum)

**Elution:** 4 x 0.5 ml CH<sub>2</sub>Cl<sub>2</sub> (percolate first 0.5 ml into the column packing without vacuum, then apply light vacuum), if necessary evaporate in a stream of nitrogen and fill up with a suitable solvent

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
### PAHs from soil

**Column type:**  
CHROMABOND® CN/SiOH,  
6 ml, 500/1000 mg

**Sample pretreatment:** dry 30 g soil with sodium sulphate and reflux 4 h with 250 ml petroleum ether in a Soxhlet extractor. For low PAH contents (colourless or weakly coloured extracts) concentrate extract to 1/10 of its volume in a rotation evaporator.

**Column conditioning:** 4 ml petroleum ether

MN Appl. Nr. 301310

**Sample application:**

aspirate 20 ml of the extract through the column


**Washing:** 2 ml petroleum ether

**Elution:** 2 x 2 ml acetonitrile / toluene (3:1, v/v), then evaporate or fill to the volume required

Further analysis: HPLC, e.g. with column 250 x 3 mm NUCLEOSIL® 5 C<sub>18</sub> PAH

recovery rates see application 301310 at [www.mn-net.com](http://www.mn-net.com)

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**Internal standard solution:** dissolve 20 mg *n*-tetracontane (C<sub>40</sub>H<sub>82</sub>) in petroleum ether, add 20 ml *n*-decane (C<sub>10</sub>H<sub>22</sub>) and fill up to 1 litre with petroleum ether. For preparation of the extraction solution dilute standard solution 1:10 with petroleum ether.

**Sample pretreatment:** adjust 900 ml water (10 °C) with HCl (12 mol/l) to pH 2 and add 80 g MgSO<sub>4</sub>. Add 50 ml of the extraction solution, close the bottle and stir the suspension intensely for 30 min.

Add enough dist. water to separate the organic from the aqueous phase.

**Column conditioning:** 5 ml petroleum ether


**Sample application:**

slowly aspirate or force the sample through the column  
**Elution:** wash with 10 ml petroleum ether. Evaporate the combined solutions from sample application and elution to 1 ml at about 75 °C. If necessary, fill up to 1 ml again. (If the hydrocarbon content is high, evaporation to 1 ml may not be necessary.)

**Recovery rate:** must be > 80 % for *n*-tetracontane.

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**Column conditioning:** 1 ml *n*-hexane

**Sample application:**

apply 250 µl waste oil to the column and aspirate or force it into the adsorbent with 2 x 1 ml *n*-hexane

MN Appl. No. 301390

**Elution:** aspirate or force another 2 x 500 µl *n*-hexane through the column; collect all *n*-hexane fractions and if necessary adjust to a concentration suitable for subsequent analysis by either evaporating *n*-hexane in a stream of nitrogen or by dilution with *n*-hexane

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
### Protection from cross contamination

For special applications, which require maximum protection from cross contamination we supply chrome-plated brass valves and stainless steel or PTFE connectors, the application of which is shown below. These special connectors are fitted through the lid; thus the sample only has contact with the inert connector and not with the lid, directly flowing into the receptacle.

### Drying attachment

If the eluate has to be evaporated, this can be performed with the so-called drying attachment (**11**, see below). This special lid has a gas connector on one side (**12**), from which the gas is fed simultaneously to the 12 or 24 stations (**13**). Thus 12 or 24 eluates can be evaporated simultaneously by just changing the lid and applying a stream of inert gas, e.g. nitrogen.