Syringe filters CHROMAFIL®

Syringe filters are used for filtration of suspended matter from liquid samples or gases. With CHROMAFIL®, rapid purification and removal of particles is very simple: just place the filter on the syringe, and you are ready for filtration. Special manipulations are not required. Contamination of sensitive instrumentation by solid impurities can be avoided, thus increasing lifetime of chromatographic columns and equipment.

Advantages:

- ◆ **Polypropylene housing**
considerably better solvent stability compared to acrylate and polystyrene filters, low content of extractable substances
- ◆ **Housing ultrasonically sealed, not glued**
no extractable components from glues
- ◆ The special **thick rim** of the housing is ideal for use of the filters in laboratory robots (e.g. Benchmate™).
- ◆ **Filtration in both directions** possible, the liquid cannot bypass the membrane
- ◆ **Luer lock on side of entry**
safe connection on the "high pressure" side
- ◆ **Luer exit**
standard luer for 25 mm filters, minispikes luer with low dead volume and small OD for 15 mm filters. Filter inlet and filter exit can be fitted to the CHROMABOND® columns for selective sample preparation with the aid of a special adaptor.
- ◆ **Deflector**
the stream of liquid is broken and distributed, and does not directly hit the membrane: this prevents rupture of the membrane
- ◆ **Star-shaped distribution device**
the liquid is evenly distributed to the whole membrane surface: this results in a better utilisation of the total area; the filter is not plugged up rapidly; high flow efficiency
- ◆ **Colour coded filters**
filters with 0.2 µm pores have a yellow upper shell, that of filters with 0.45 µm pores is colourless; the different membrane types are distinguished by different colours of the lower shell
- ◆ Available **pore sizes** 0.2 and 0.45 µm (exceptions: PET filters with 1.2 µm, glass fibre filters with 1 µm, PES filters with 5 µm)
- ◆ **Filter sizes:** 25 and 15 mm diameter The small diameter filters are especially recommended for very small samples, which require extremely low dead volumes: 80 µl for 25 mm Ø, 12 µl for 15 mm Ø
- ◆ All filters can be **autoclaved** at **121 °C** and **1.1 bar** for 30 min.

Recommended filter size depending on sample volume

| sample volume | recommended filter diameter |
|---------------|-----------------------------|
| 1 - 10 ml | 15 mm |
| 10 - 100 ml | 25 mm |

CHROMAFIL® BIG-BOXES

- ◆ 400 (25 mm) or 800 (15 mm) quality syringe filters
- ◆ food safe PE box with screw cap
- ◆ economical price

Depending on your filtration task you can choose filter membranes made from different materials:

| Material |
|----------------------------------------------------------------------------------|
| Polyester (PET) with or without glass fibre prefilter |
| Regenerated cellulose (RC) |
| Teflon® (PTFE) |
| Cellulose mixed esters (MV) |
| Cellulose acetate (CA) · sterile and non-sterile |
| Polyamide / Nylon (PA) |
| Polyethersulfone (PES) · sterile |
| Polyvinylidene difluoride (PVDF) with or without glass fibre prefilter |
| Glass fibre (GF) |

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Chemical compatibility of filter materials

The following table lists the chemical compatibility of our CHROMAFIL® materials. The chemical compatibility depends on several parameters such as time, pressure, temperature and concentration.

In most cases, CHROMAFIL® filters will have only short contact with a solvent. In these cases they may be used despite of limited compatibility.

For example, a PTFE filter with PP housing does not liberate any UV-detectable substances during filtration of 5 ml THF, although PP shows only limited resistance towards THF.

| Solvent | Material | | | | | | | | | |
|------------------------------|----------|----|----|----|------|------|-----|-----|----|----|
| | MV | CA | RC | PA | PTFE | PVDF | PES | PET | GF | PP |
| Acetaldehyde | ⊖ | ⊖ | ⊕ | ⊙ | ⊕ | ⊕ | | ⊕ | ⊕ | ⊙ |
| Acetic acid, 100 % | ⊖ | ⊖ | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Acetone | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ |
| Acetonitrile | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Ammonia, 25 % | ⊖ | ⊖ | ⊙ | ⊖ | ⊕ | ⊕ | ⊕ | ⊙ | ⊕ | ⊕ |
| Benzene | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊙ | | ⊕ | ⊕ | ⊙ |
| n-Butanol | ⊕ | ⊕ | ⊕ | ⊙ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Cyclohexane | ⊕ | ⊕ | ⊕ | ⊙ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Dichloromethane | ⊕ | ⊖ | ⊕ | ⊖ | ⊕ | ⊕ | ⊖ | ⊕ | ⊕ | ⊖ |
| Diethyl ether | ⊙ | ⊙ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊙ |
| Dimethylformamide | ⊖ | ⊖ | ⊙ | ⊕ | ⊕ | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ |
| 1,4-Dioxane | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ | ⊙ | ⊖ | ⊕ | ⊕ | ⊙ |
| Ethanol | ⊖ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Ethyl acetate | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊙ |
| Ethylene glycol | ⊙ | ⊙ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Formic acid, 100 % | ⊕ | ⊖ | ⊙ | ⊖ | ⊕ | ⊕ | ⊕ | ⊙ | ⊕ | ⊕ |
| Hydrochloric acid, 30 % | ⊖ | ⊖ | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ | ⊖ | ⊕ | ⊕ |
| Methanol | ⊖ | ⊖ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Nitric acid, 65 % | ⊖ | ⊖ | ⊖ | ⊖ | ⊙ | ⊙ | | ⊙ | ⊕ | ⊖ |
| Oxalic acid, 10 % aqueous | ⊕ | ⊖ | ⊕ | ⊖ | ⊕ | ⊕ | | ⊕ | ⊕ | ⊕ |
| Petroleum ether | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Phosphoric acid, 80 % | ⊖ | ⊖ | ⊙ | ⊖ | ⊕ | ⊙ | | ⊕ | ⊕ | ⊕ |
| Potassium hydroxide, 1 mol/l | ⊖ | ⊖ | ⊙ | ⊕ | ⊕ | ⊙ | ⊕ | ⊙ | ⊕ | ⊕ |
| 2-Propanol | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Sodium hydroxide, 1 mol/l | ⊖ | ⊖ | ⊙ | ⊕ | ⊕ | ⊙ | ⊙ | ⊙ | ⊙ | ⊕ |
| Tetrachloromethane | ⊕ | ⊖ | ⊕ | ⊕ | ⊕ | ⊙ | | ⊕ | ⊕ | ⊙ |
| Tetrahydrofuran | ⊖ | ⊖ | ⊕ | ⊙ | ⊕ | ⊕ | ⊖ | ⊕ | ⊕ | ⊙ |
| Toluene | ⊕ | ⊖ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊙ |
| Trichloroethene | ⊕ | ⊕ | ⊕ | ⊙ | ⊕ | ⊕ | | ⊕ | ⊕ | ⊙ |
| Trichloromethane | ⊕ | ⊖ | ⊕ | ⊖ | ⊕ | ⊕ | ⊖ | ⊕ | ⊕ | ⊖ |
| Urea | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | | ⊕ | ⊕ | ⊕ |
| Water | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ |
| Xylene | ⊕ | ⊕ | ⊕ | ⊕ | ⊕ | ⊙ | | ⊕ | ⊕ | ⊙ |

Data not guaranteed. ⊕ resistant, ⊖ not resistant, ⊙ limited resistance

MV = cellulose mixed esters, CA = cellulose acetate, RC = regenerated cellulose, PA = polyamide, PTFE = polytetrafluoroethylene (Teflon), PVDF = polyvinylidene difluoride, PES = polyethersulfone, PET = polyester, GF = glass fibre, PP = polypropylene (housing material)

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

Xtra: polipropilene piano colour-free.

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